

USER MANUAL

Diabetes Risk Analyzer v1.0

USER MANUAL

Diabetes Risk

TABLE OF CONTENTS

1. Introduction

2. System Re

1. INTRODUCTION

Diabetes Risk Analyzer is an Android application designed to

estimate your

The app analyzes the following inputs:

- Patient nam

2. SYSTEM REQUIREMENTS

- Android Version: 7.0 (Nougat) or higher

- RAM: Minim

3. INSTALLATION GUIDE

METHOD A ? Install from APK file:

Step 1: Download the APK file from the GitHub release page

Step 2: Open

METHOD B ? Run from Android Studio:

Step 1: Clone the project from GitHub

Step 2: Open

4. APP SCREENS OVERVIEW

The app contains the following screens:

SCREEN 1 ? Splash Screen

Shows the a

SCREEN 2 ? Main Input Screen (Home)

Where you

SCREEN 3 ? Result Screen

Shows your

SCREEN 4 ? History Screen

Shows all y

5. STEP-BY-STEP USAGE GUIDE

STEP 1 ? OPEN THE APP

Launch Dial

STEP 2 ? ENTER PATIENT INFORMATION

Full Name:

STEP 3 ? ANSWER 5 SYMPTOM QUESTIONS

Q1 - Daily Water Intake:

For each qu

Q2 - Urination Frequency:

Select h

Q3 - Weight Loss (Last 6 Months):

Select h

Q4 - Wound Healing Time:

Select h

Q5 - Hunger Level:

Select h

NOTE: The progress bar at the top shows how many

Select h

questions y

STEP 4 ? TAP ANALYZE

If any field is empty, an error message will appear

Once all fiel

pointing to t

STEP 5 ? VIEW YOUR RESULTS

SUMMARY TAB shows:

The Result

- Risk level

CHART TAB shows:

- Bar chart c

STEP 6 ? SHARE YOUR RESULT (Optional)

Tap the Sha

STEP 7 ? VIEW HISTORY

Tap the His

6. UNDERSTANDING YOUR RESULTS

RISK LEVELS EXPLAINED:

GREEN ? LOW RISK (Score 7 or below)

Your respon

YELLOW ? MEDIUM RISK (Score 8 to 14)

Some warni

RED ? HIGH RISK (Score 15 or above)

Multiple sign

SCORE BREAKDOWN:

Base Score

7. HISTORY FEATURE

The app automatically saves every analysis you complete.

To VIEW history:

Tap the His

Each history record shows:

- Patient na

To DELETE history:

Open Histor

NOTE: History is stored locally on your device.

Uninstalling th

8. NOTIFICATIONS

After every analysis, the app sends a push notification

to your device

Notification types:

LOW RISK:

To enable notifications (Android 13+):

When prom

9. TROUBLESHOOTING

PROBLEM: App crashes on opening

SOLUTION: F

PROBLEM: Analyze button not working

SOLUTION: M

PROBLEM: No notification appearing

SOLUTION: C

PROBLEM: History screen is empty

SOLUTION: Y

PROBLEM: Charts not showing

SOLUTION: S

PROBLEM: BMI field ? what value to enter

SOLUTION: C

10. MEDICAL DISCLAIMER

IMPORTANT ? PLEASE READ CAREFULLY

This application is developed for EDUCATIONAL and AWARENESS

purposes ON

The Diabetes Risk Analyzer:

- Does NOT p

The results provided by this app are ESTIMATES only and

may not accu

(doctor, endo

extreme thirst

For support o

ALWAYS consult a qualified, licensed medical professional

If you are experiencing symptoms of diabetes such as

Thank you for using Diabetes Risk Analyzer.